



product details

QTP-M 2X26-32

Product description: QTP-M 2X26-32
Product code: 4008321329165
Quantity: Shipping carton box (VS) contains 20 Piece (PCE)

You can find this product in the eCatalog:

http://catalog.myosram.com?~language=EN&~country=DE&it_p=4008321329165



Applications	
Suitable f. fixtures with prot. class	Class I Class II
Categorizations	
SEG number	7980317
Certificates & Standards	
Approval marks - approval	ENEC 10 VDE VDE-EMC
EEL - Energy Label	A2
Standards	acc. to EN 61347-2-3 acc. to EN 60929 acc. to EN 55015:2006 + A1:2007 acc. to EN 61000-3-2 acc. to EN 61547
Type of protection	IP20
General Description	
Design / version	ECG multiwatt
New article	NEW
Logistical Data	
Product weight	180 g
Relative humidity	5...85 % ¹⁾
Temperature range at storage	-40...85 °C
Technical - Electrical Data	
Maximum mains current	0.30 A
Nominal voltage	220-240 V
Mains frequency	50-60 Hz
DC voltage	176...276 V
Operating frequency	45 kHz ²⁾
Output voltage	430 V
Dielectric strength	300 V 320 V 350 V ³⁾
Max. no. of ECGs on circuit breaker 10 A	17
Max. no. of ECGs on circuit breaker 16 A	28
TH full width at half maximum	250 µs
Technical - Geometries	
Length	123.00 mm
Width	79.00 mm
Nominal voltage	33.00 mm
Mounting hole spacing, length	129.5 mm
Mounting hole spacing, width	67 mm
Technical - Light Technical Data	
Rated starting time	preheat start within 1.0 s



Technical - Temperatures	
Ambient temperature range	-20...+50 °C
Operating temperature	70 °C
Max. housing temperature in case of fault	110 °C

Packaging units				
Product code	Packaging type and content	Dimensions in h x w x l	Gross weight	Volume
4008321329158	Unpacked contains 1 Piece	80,000 mm x 33,000 mm x 123,000 mm	180,000 g (0,000 g)	0,327 Cubic dec.
4008321329165	Shipping carton box contains 20 Piece	257,000 mm x 88,000 mm x 428,000 mm	3.860,000 g (0,000 g)	8,580 Cubic dec.

- 1) Maximum 56 days/year at 85 %
- 2) Approximately
- 3) Permanent | 48 h | 2 h

**Lamp wiring:**

³⁾ Max. permitted cable length in between ECG and lamp: 2.0 m (PIN 21; 22); 1.0 m (PIN 23; 24)

⁴⁾ Max. permitted cable length in between ECG and lamp: 2.0 m (PIN 21; 22); 1.0 m (PIN 25; 26)